

**1 Identification of the substances/ mixture and of the company/ undertaking****1.1 Product Identifiers**

Product Number CR027  
Product Name Yeast Extract Powder, Certified  
REACH Registration Number This product is a mixture. Reach registration number is not available for this mixture.

**1.2 Relevant identified uses of the substance or mixture and uses advised against**

**1.2.1** Relevant identified uses Laboratory Chemicals, Analytical Purpose, Biochemical Analysis

**1.3 Details of the supplier of the safety data sheet**

Produced by HiMedia Laboratories Private Limited  
Address 23, Vadhani Industrial Estate, LBS Marg, Ghatkopar (W), Mumbai - 400 086 India  
Tel. No. +91-22-2500 0970, +91-22-2500 1607 Fax No. : +91-22-25002468  
Mail Id [info@himedialabs.com](mailto:info@himedialabs.com) Website : [www.himedialabs.com](http://www.himedialabs.com)

**1.4 Emergency Tel. No.**

Emergency Tel. No. Please contact the regional HiMedia representation in your country

**2 Hazards Identification****2.1 Classification of the substance or mixture**

***CLP Classification-Regulation (EC) No. 1272/2008[EU-GHS/CLP]***

Not a hazardous substance or mixture according to Regulation (EC) No.1272/2008.

**2.2 Label elements**

***Labeling according to Regulation (EC) No.1272/2008***

The product does not need to be labelled in accordance with EC directives or respective national laws.

**2.3 Other Hazards**

None

**3 Composition/Information On Ingredients****3.2 Mixture**

The components of this mixture need not be disclosed as per the regulations.

**4 First Aid Measures****4.1 Description of first aid measures**

***General advice***

Consult a physician. Show this safety data sheet to the doctor in attendance.

***If inhaled***

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

***In case of skin contact***

Wash off with soap and plenty of water. Consult a physician.

***In case of eye contact***

Rinse immediately with plenty of water for at least 15 minutes. Consult a physician.

***If swallowed***

Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician if feeling unwell.

**4.2 Most important symptoms and effects, both acute and delayed**

No data available.

**4.3 Indication of immediate medical attention and special treatment needed**

No data available

---

**5 Fire Fighting Measures**

**5.1 Extinguishing media**

***Suitable extinguishing media***

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

***Unsuitable extinguishing media***

No data available.

**5.2 Special hazards arising from the substance or mixture**

Nature of decomposition products unknown

**5.3 Precautions for fire-fighters**

Wear full protective clothing and self contained breathing apparatus for fire fighting if necessary.

**5.4 Further information**

No data available

---

**6 Accidental Release Measures**

**6.1 Personal precautions, protective equipment and emergency procedures**

Avoid breathing vapours, mist or gas. Wear respiratory protection. Ensure adequate ventilation. Evacuate personnel to safe areas.

**6.2 Environmental precautions**

Prevent further leakage or spillage if safe to do so. Do not let product enter drains.

**6.3 Methods and materials for containment and cleaning up**

Soak up with inert absorbent material. Place material in a suitable, closed containers for disposal. Dispose off properly

**6.4 Reference to other sections**

For disposal see Section 13.

---

**7 Handling and Storage**

**7.1 Precautions for safe handling**

Avoid contact with skin and eyes. Avoid inhalation of vapour or mist. Normal measures for preventive fire protection.

**7.2 Conditions for safe storage, including any incompatibilities**

Store in cool place. Keep container tightly closed in a dry and well-ventilated place. Containers which are opened must be carefully resealed and kept upright to prevent leakage.

**Recommended Storage Temperature** : On receipt store between 10-30°C

### 7.3 Specific end uses

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated.

---

## 8 Exposure Controls/Personal Protection

### 8.1 Control parameters

Components with workplace control parameters

### 8.2 Exposure controls

#### **Appropriate engineering controls**

Avoid contact with skin, eyes and clothing. Wash hands before breaks and immediately after handling the products.

#### **Personal protective equipment**

##### **Hygiene measure**

Immediately change contaminated clothing. Apply preventive skin protection. Wash hands and face after working with the product.

##### **Eye/face protection**

Faceshield and safety goggles. Faceshield (8-inch minimum). Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166 (EU). Have eye-washing facilities readily available where eye contact can occur.

##### **Skin protection**

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands. The selected protective gloves have to satisfy the specifications of EU Directive 89/686/EEC and the standard EN 374 derived from it.

##### **Body protection**

Complete suit protecting against chemicals. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

##### **Respiratory protection**

Where risk assessment shows air-purifying respirators are appropriate use a full-face respirator with multi-purpose combination (US) or type ABEK (EN 14387) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

##### **Environment exposure controls**

Do not empty into drains.

---

## 9 Physical and chemical properties

### 9.1 Information on basic physical and chemical properties

|                        |   |
|------------------------|---|
| Appearance             | Light yellow to brownish yellow coloured homogenous free flowing powder |
| Odour                  | Characteristic odour but not putrescent                                 |
| Odour Threshold        | No data available   |
| pH                     | 6.50 - 7.50   |
| Melting/freezing point | No data available   |

|   |  |
|---|--|
| Initial boiling point and boiling range | No data available  |
| Flash point                             | No data available  |
| Flammability (Solid, gas)               | No data available  |
| Vapour pressure                         | No data available  |
| Relative density                        | No data available  |
| Water Solubility                        | Soluble in distilled/ purified water, insoluble in alcohol and ether |
| Partition coefficient: n-octanol/water  | No data available  |
| Autoignition Temperature                | No data available  |
| Decomposition Temperature               | No data available  |
| Explosive properties                    | No data available  |
| Oxidizing properties                    | No data available  |
| Vapour density                          | No data available  |
| Thermal decomposition                   | No data available  |

## 9.2 Other safety information

No data available

## 10 Stability and Reactivity

### 10.1 Reactivity

No data available

### 10.2 Chemical stability

No data available

### 10.3 Possibility of hazardous reactions

No data available

### 10.4 Conditions to avoid

No data available

### 10.5 Incompatible materials

No data available

### 10.6 Hazardous decomposition products

Refer Section 5.2

## 11 Toxicological Information

### 11.1 Information on toxicological effects

#### *Acute toxicity*

No data available

#### *Skin corrosion/irritation*

No data available

#### *Serious eye damage/eye irritation*

No data available

#### *Respiratory or skin sensitisation*

No data available

#### *Germ cell mutagenicity*

No data available

#### *Carcinogenicity*

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

**Reproductive toxicity**

No data available

**Specific target organ toxicity- single exposure**

No data available

**Aspiration hazard**

No data available

**Potential Health Effects**

**Inhalation**

REFER SECTION 2

**Skin**

REFER SECTION 2

**Eyes**

REFER SECTION 2

**Ingestion**

REFER SECTION 2

**Additional Information**

RTECS : No Data Available

---

**12 Ecological Information**

**12.1 Toxicity**

No data available

**12.2 Persistence and degradability**

No data available

**12.3 Bioaccumulative potential**

No data available

**12.4 Mobility in soil**

No data available

**12.5 PBT and vPvB assessment**

This substance or mixture contains no components considered to be persistent, bioaccumulating nor toxic (PBT) at levels of 0.1% or higher.

**12.6 Other adverse effects**

No data available

---

**13 Disposal Considerations**

**13.1 Waste treatments methods**

**Product**

Offer surplus and non-recyclable solutions to a licenced disposal company. Contact a licenced professional waste disposal service to dispose off this material.

**13.2 Contaminated packaging**

Dispose of as unused product.

## 14 Transport Information

### 14.1 UN-No

ADNR : ADR : IATA\_C : IATA\_P : IMDG : RID :

### 14.2 UN proper shipping name

ADNR : Not dangerous goods  
ADR : Not dangerous goods  
IATA\_C : Not dangerous goods  
IATA\_P : Not dangerous goods  
IMDG : Not dangerous goods  
RID : Not dangerous goods

### 14.3 Transport hazard class(es)

ADNR : - ADR : - IATA\_C : - IATA\_P : - IMDG : - RID : -

### 14.4 Packaging group

ADNR : ADR : IATA\_C : IATA\_P : IMDG : RID :

### 14.5 Environmental hazards

ADNR : NO ADR : NO IMDG : Marine pollutant - NO IATA\_C : NO IATA\_P : NO RID : NO

### 14.6 Special precautions for use

No data available

## 15 Regulatory Information

This safety data sheet complies with the requirements of Regulation (EC) No. 1907/2006

### 15.1 Safety health and environment regulations/legislation specific for the substance or mixture

No data available

### 15.2 Chemical Safety Assessment

No data available

## 16 Other information

### Further Information

Copyright 2016 HiMedia Laboratories Pvt. Ltd.

The information given in this safety data sheet is believed to be correct yet does not claim to be all inclusive. This document is intended only as a guide for appropriate precautionary handling of the material by properly trained individuals, information here being commensurate with the present state of our knowledge regarding the manner and conditions of use, handling, storage or disposal. The information provided herein shall not be considered as guarantee of the properties of the product. HiMedia Laboratories, shall not be held liable for any damage resulting from improper handling or contact with the above product. Unless explicitly stated on the product or in any of the documentation accompanying the product, it is intended for research or testing purpose only and is not to be used for any other purpose.

