
Sodium acetate anhydrous, Hi-AR™

GRM1012

Product Identifier

CAS No.	:	127-09-3
EC No.	:	204-823-8
Molecular Formula	:	C ₂ H ₃ O ₂ Na
Molecular Weight	:	82.03
Synonym	:	Acetic acid sodium salt
HS Code	:	2915 29 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white hygroscopic crystals or powder or solid
Solubility	:	100 mg soluble in 1 mL of water
pH (5% in water at 25°C)	:	7.00 - 9.20
Chloride (Cl)	:	<= 0.002%
Heavy metals (as Pb)	:	<= 0.001%
Water insoluble matter	:	<= 0.01%
Iron (Fe)	:	<= 0.001%
Phosphate (PO ₄)	:	<= 0.001%
Sulfate (SO ₄)	:	<= 0.003%
Loss on drying (at 120°C, 2 hr)	:	<= 1.0%
Assay (NT, on dry basis)	:	min. 99.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	AJ4300010
WGK	:	1