

---

## Magnesium sulphate dried, Hi-AR™

**GRM1281**

### Product Identifier

---

CAS No.	:	7487-88-9
EC No.	:	231-298-2
Molecular Formula	:	MgSO <sub>4</sub>
Molecular Weight	:	120.37
HS Code	:	2833 21 00
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	White to light grey crystals or granules or powder
Solubility	:	33.3 mg soluble in 1 mL of hot water
pH (5% in water at 25°C)	:	5.00 - 9.20
Loss on ignition (at 105°C, 2hr then ignite at 450 ±25°C, 2hr)	:	22.00 - 28.00%
Assay (EDTA Titration, on ignited basis)	:	99.00 - 102.00%

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	OM4500000
WGK	:	1