

GRM1311 Potassium sorbate**Product Number Packing**

GRM1311	: 1KG
GRM1311	: 500G

Product Information

Product Code	: GRM1311
Product Name	: Potassium sorbate
Synonym	: Sorbic acid potassium salt; Potassium 2,4-hexadienoate
Molecular Formula	: C ₆ H ₇ KO ₂
Molecular Weight	: 150.22
CAS No.	: 24634-61-5
EC No.	: 246-376-1
HS Code	: 2916 19 95
Shelf Life	: 4 years

Technical Specification


Appearance	: White to light cream granules or powder or beads
Solubility	: 1000 mg soluble in 1 mL of water
Alkalinity (as K ₂ CO ₃)	: <= 1.00%
Acidity	: <= 1.00%
Loss on drying (at 105°C, 3hr)	: <= 1.00%
Assay (NT, on dry basis)	: 99.00 - 101.00 %

GHS Safety Information

Hazard Statement(s)	: H315-H319-H335
Precautionary Statement(s)	: P261-P305+P351+P338
Signal Word	: Warning
Hazard Pictogram(s)	

**Risk and Safety Information**

R-Phrase(s)	: 36/37/38
S-Phrase(s)	: 26
WGK	: 1
RTECS	: WG2170000
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	

**Transport Information**

Marine Pollutant	: No
------------------	------



Product Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein