

## Material Safety Data Sheet

PAGE 1 OF 2

Name of the Product : L (+) Tartaric acid, Extra pure

Code No. : **RM1430** 

Section 1 : Chemical Identification

Code No. : RM1430

Name of the Product : L (+) Tartaric acid, Extra pure Produced by : HiMedia Laboratories Pvt. Ltd.

Address : 23, Vadhani Indl. Estate, LBS Marg, Mumbai 400 086, India.

S

Tel. No. : 2500 0970, 2500 1607 Fax No. 022 2500 2468

Section 2 : L(+) Tartaric acid, Extra pure

Irritant material

L (+) Tartaric acid, Extra pure

CAS No. : 87-69-4 R : 36/37/38

RTECS : WW 7875000

Section 3 : Hazards Identification

Hazard: This symbol designates substances which may have an irritant effect on skin, eyes and respiratory

: 26-36

organs.

Caution: Do not breathe vapours and avoid contact with skin and eyes.

Section 4 : First - Aid Measures

In case of swallowing, wash out mouth with water provided person is conscious. In case of skin contact, flush with copious amounts of water for at least 15 minutes. Remove contaminated clothing and shoes. In case of inhalation, remove to fresh air. If breathing becomes difficult, call physician. In case of contact with eyes, wash with copious amounts of water for at least 15 minutes. Assure adequate washing by separating the eyelids with fingers. In serious conditions, call a physician, show the container or label.

Section 5 : Fire Fighting Measures

Extinguishing Media. Water spray.

Carbon dioxide, dry chemical powder or appropriate foam. Special firefighting procedures. Wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.

Section 6 : Accidental Release Measures

Wear respirator, chemical safety goggles, rubber boots and heavy rubber gloves. Sweep up, place in a bag and hold for waste disposal. Avoid raising dust. Ventilate area and wash spill site after material pickup is complete.

Section 7 : Handling and Storage

Handling - Refer to Section 8 Storage - Store below 30°C

Section 8 : Exposure Controls/Personal Protection

Wear appropriate NIOSH/MSHA-approved respiratory, chemical-resistant gloves, safety goggles, other protective clothing. Mechanical exhaust required.



## Material Safety Data Sheet

PAGE 2 OF 2

Section 9 : Physical and Chemical Properties

Appearance : Crystalline powder

Colour : White

Section 10 : Stability and Reactivity

Stability: Product is stable if stored as per the conditions specified under storage of Section No. 7. Product

loses its potency/performance above 45°C.

Conditions to avoid: Heat and light. Hazardous polymerization will not occur.

Section 11 : **Toxicological Information** 

Acute Effects: May be harmful if swallowed.

Exposure can cause: Stomach pains, vomiting, diarrhoea. Prolonged or repeated exposure may

cause allergic reactions in certain sensitive individuals.

RTECS No.: WW 7875000

For complete information RTECS may be referred

Section 12 : Ecological Information

Data not available

Section 13 : Disposal Considerations

On completion of work all used or unusable preparations of this product and derivatives thereof are to be disposed off by autoclaving and/or by incineration. Dispose off waste in accordance with all applicable

Federal, State and local laws.

Section 14 : Transport Information

UN No.: Not applicable.

Section 15 : Regulatory Information

Risk Phrases : Irritating to eyes, respiratory system and skin

Safety Phrases: In case of contact with eyes, rinse immediately with plenty of water and seek medical

advice

Safety Phrases: Wear suitable protective clothing

Significance of signs: This symbol designates substances which may have an irritant effect on skin, eyes

and respiratory organs.

Section 16 : Other Information

The information contained in this data sheet represents the best information currently available to us. However, no warranty is made with respect to its completeness and we assume no liability resulting from its use. The information is offered solely for users' obligation to investigate and determine the suitability of the

information for their particular purpose.