

1 Identification of the substances/ mixture and of the company/ undertaking**1.1 Product Identifiers**

Product Number	GRM1880
Product Name	2,4-Hexadienoic acid, Hi-LR™
CAS No. :	110-44-1
EC No. :	203-768-7

1.2 Relevant identified uses of the substance or mixture and uses advised against**1.2.1** Relevant identified uses Laboratory chemicals, Manufacture of substances**1.3 Details of the supplier of the safety data sheet**

Produced by	HiMedia Laboratories Private Limited		
Address	23, Vadhani Industrial Estate, LBS Marg, Ghatkopar (W), Mumbai - 400 086 India		
Tel. No.	+91-22-2500 0970, +91-22-2500 1607	Fax No. :	+91-22-25002468
Mail Id	info@himedialabs.com	Website :	www.himedialabs.com

1.4 Emergency Tel. No.

Emergency Tel. No. Please contact the regional HiMedia representation in your country

2 Hazards Identification**2.1 Classification of the substance or mixture****CLP Classification-Regulation (EC) No. 1272/2008[EU-GHS/CLP]**

Skin corrosion or irritation, (Category 2), H315

Serious eye damage or eye irritation, (Category 2A), H319

Specific target organ toxicity, single exposure, Respiratory tract irritation, (Category 3), H335

2.2 Label elements**Labeling according to Regulation (EC) No.1272/2008**Pictogram
Signal word Warning

Hazard Statement(s)

H315	Causes skin irritation
H319	Causes serious eye irritation
H335	May cause respiratory irritation

Precautionary Statement(s)

P261	Avoid breathing dust/fume/gas/mist/vapours/spray.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 IF eye irritation persists: Get medical advice/attention.

2.3 Other Hazards

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

3 Composition/Information On Ingredients

3.1 Substances

Synonym	Sorbic acid, Hi-LR™
Molecular Formula :	C ₆ H ₈ O ₂
Molecular Weight :	112.13
CAS No. :	110-44-1
EC No. :	203-768-7

Component	Classification	Concentration
2,4-Hexadienoic acid, Hi-LR™		
CAS No. : 110-44-1 EC No. : 203-768-7	As Per EC Regulation 1272/2008 Skin Irrit. 2; Eye Irrit. 2A; STOT SE 3 H315; H319; H335	100%

For the full text of the H-Statements and classification mentioned in this Section, see Section 16

4 First Aid Measures

4.1 Description of first aid measures

General advice

Consult a physician. Show this safety data sheet to the doctor in attendance.

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

In case of skin contact

Wash with plenty of soap and water. Consult a physician.

In case of eye contact

Rinse immediately with plenty of water for at least 15 minutes. Consult a physician.

If swallowed

Do NOT induce vomiting. Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

4.2 Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in the labeling (see section 2.2) and/or in section 11.

4.3 Indication of immediate medical attention and special treatment needed

No data available

5 Fire Fighting Measures

5.1 Extinguishing media

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Unsuitable extinguishing media

No data available.

5.2 Special hazards arising from the substance or mixture

Carbon oxides

5.3 Precautions for fire-fighters

Wear self contained breathing apparatus for fire fighting if necessary

5.4 Further information

No data available

6 Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear respiratory protection. Avoid breathing vapours, mist or gas. Ensure adequate ventilation.

Evacuate personnel to safe areas.

6.2 Environmental precautions

Prevent further leakage or spillage if safe to do so. Do not let product enter drains. Discharge into environment must be avoided.

6.3 Methods and materials for containment and cleaning up

Sweep up and shovel. Contain spillage, and then collect with an electrically protected vacuum cleaner or by wet-brushing and place in container for disposal according to local regulations (see section 13). Keep in suitable, closed containers for disposal.

6.4 Reference to other sections

For disposal see Section 13.

7 Handling and Storage

7.1 Precautions for safe handling

Avoid contact with skin and eyes. Avoid inhalation of dust and aerosols. Provide appropriate exhaust ventilation at places where dust is formed. For precautions see section 2.2

7.2 Conditions for safe storage, including any incompatibilities

Store in cool place. Keep container tightly closed in a dry and well-ventilated place.

Recommended Storage Temperature : Below 30°C

7.3 Specific end uses

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated.

8 Exposure Controls/Personal Protection

8.1 Control parameters

Components with workplace control parameters

8.2 Exposure controls

Appropriate engineering controls

Handle in accordance to general industrial hygiene and safety practice. Wash hands before breaks and at the end of workday.

Personal protective equipment

Hygiene measure

Immediately change contaminated clothing. Apply preventive skin protection. Wash hands and face after working with the product.

Eye/face protection

Faceshield and safety goggles.

Skin protection

Handle with gloves.

Body protection

Complete suit protecting against chemicals.

Respiratory protection

Air-purifying respirators.

Environment exposure controls

Prevent further leakage or spillage if safe to do so. Do not let product enter drains. Discharge into the environment must be avoided.

9 Physical and chemical properties

9.1 Information on basic physical and chemical properties

Appearance	White to cream crystals or powder
Odour	No data available
Odour Threshold	No data available
pH	No data available
Melting/freezing point	No data available
Initial boiling point and boiling range	No data available
Flash point	No data available
Upper/lower flammability or explosive limits	No data available
Evaporation rate	No data available
Flammability (Solid, gas)	No data available
Vapour pressure	No data available
Relative density	No data available
Water Solubility	No data available
Partition coefficient: n-octanol/water	No data available
Autoignition Temperature	No data available
Decomposition Temperature	No data available
Viscosity	No data available
Explosive properties	No data available
Oxidizing properties	No data available
Vapour density	No data available
Thermal decomposition	No data available

9.2 Other safety information

No data available

10 Stability and Reactivity

10.1 Reactivity

- No data available
- 10.2 Chemical stability**
Stable under recommended storage conditions.
- 10.3 Possibility of hazardous reactions**
No data available
- 10.4 Conditions to avoid**
Incompatible material
- 10.5 Incompatible materials**
Bases, Oxidizing agents, Strong oxidizing agents, Strong reducing agents
- 10.6 Hazardous decomposition products**
Other decomposition products - refer section 5.2.
In the event of fire - refer section 5.

11 Toxicological Information

11.1 Information on toxicological effects

Acute toxicity

No data available

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Carcinogenicity

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

Reproductive toxicity

No data available

Specific target organ toxicity- single exposure

No data available

Specific target organ toxicity - repeated exposure

No data available

Aspiration hazard

No data available

Additional Information

RTECS : WG2100000

12 Ecological Information

12.1 Toxicity

No data available

12.2 Persistence and degradability

No data available

12.3 Bioaccumulative potential

No data available

12.4 Mobility in soil

No data available

12.5 PBT and vPvB assessment

PBT and vPvB assessment not available Chemical safety assessment is not required

12.6 Other adverse effects

Discharge into the environment must be avoided.

13 Disposal Considerations**13.1 Waste treatments methods****Product**

Offer surplus and non-recyclable solutions to a licensed disposal company. Dissolve or mix the material with a combustive solvent and burn in chemical incinerator equipped with an afterburner and scrubber.

13.2 Contaminated packaging

Dispose of as unused product.

14 Transport Information**14.1 UN-No**

ADNR : ADR : IATA_C : IATA_P : IMDG : RID :

14.2 UN proper shipping name

ADNR : Not dangerous goods
ADR : Not dangerous goods
IATA_C : Not dangerous goods
IATA_P : Not dangerous goods
IMDG : Not dangerous goods
RID : Not dangerous goods

14.3 Transport hazard class(es)

ADNR : - ADR : - IATA_C : - IATA_P : - IMDG : - RID : -

14.4 Packaging group

ADNR :- ADR :- IATA_C :- IATA_P :- IMDG :- RID :-

14.5 Environmental hazards

ADNR : No ADR : No IMDG : Marine Pollutant : No IATA_C : No IATA_P : No RID : No

14.6 Special precautions for use

No data available

15 Regulatory Information

This safety datasheet complies with the requirements of Regulation(EC) No. 1907/2006.

15.1 Safety health and environment regulations/legislation specific for the substance or

mixture

No data available

15.2 Chemical Safety Assessment

No data available

16 Other information

Text of H-Statements and EUH-Phrases mentioned in section 3

H315	Causes skin irritation
H319	Causes serious eye irritation
H335	May cause respiratory irritation
Skin Irrit. 2	Skin corrosion or irritation, Category 2
Eye Irrit. 2A	Serious eye damage or eye irritation, Category 2A
STOT SE 3	Specific target organ toxicity, single exposure, Respiratory tract irritation, Category 3

Further Information

Copyright 2016 HiMedia Laboratories Pvt. Ltd.

The information given in this safety data sheet is believed to be correct yet does not claim to be all inclusive. This document is intended only as a guide for appropriate precautionary handling of the material by properly trained individuals, information here being commensurate with the present state of our knowledge regarding the manner and conditions of use, handling, storage or disposal. The information provided herein shall not be considered as guarantee of the properties of the product. HiMedia Laboratories, shall not be held liable for any damage resulting from improper handling or contact with the above product. Unless explicitly stated on the product or in any of the documentation accompanying the product, it is intended for research or testing purpose only and is not to be used for any other purpose.