

---

## M-Protein, Vitamin Free

**GRM497**

---

### Product Identifier

CAS No.	:	9000-71-9
EC No.	:	232-555-1
HS Code	:	3501 10 00
Storage	:	Below 30°C
Shelf life	:	4 years

---

### Technical Specification

Appearance	:	White to cream granules or powder
Solubility	:	33.3 mg soluble in 1 mL of 1M sodium hydroxide forming hazy/cloudy solution
Vitamin content	:	Passes (determined by Microbiological methods as per AOAC)
Residue on ignition	:	<= 2.50%
Loss on drying (at 105°C, 3 hr)	:	<= 12.00%
Fat content	:	<= 1.50%
Total nitrogen (On dry basis)	:	13.00 - 16.00%
Protein content (Based on total nitrogen, dry basis)	:	min. 85.00%

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3