

GRM7494 Sodium alginate

Product Number Packing

GRM7494 : 500G
GRM7494 : 5KG

Product Information

Product Code : GRM7494
Product Name : Sodium alginate
Synonym : Alginic acid sodium salt
CAS No. : 9005-38-3
HS Code : 3913 10 00
Other Information : Not suitable for bead formation
Shelf Life : 4 years

Technical Specification

Appearance : White to yellowish fibrous or granular powder
Solubility : Soluble in water forming a viscous, colloidal solution
Identification test for alginate : Passes test
Loss on drying (at 105°C, 4 hr) : <= 15.00%

Risk and Safety Information

WGK : 1
Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods