

---

## Casein, (Hammersten)

**RM10675**

---

### Product Identifier

CAS No.	:	9000-71-9
EC No.	:	232-555-1
HS Code	:	3501 10 00
Storage	:	Below 30°C
Shelf life	:	4 years

---

### Technical Specification

Appearance	:	White to cream powder or granules
Solubility	:	33.3 mg soluble in 1 mL of alkali hydroxide may give hazy solution
Sulphated ash (800°C)	:	<= 1.50%
Loss on drying (at 105°C, 1 hr)	:	<= 12.00%
Assay (Nitrogen, on dry basis)	:	min. 95.00%

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available